

YOUR PATH TO PRISTINE PROTECTION

EMC4 GLE-582

AUTOCLAVE



Best Quality



Best Price




Warranty



Sterilization and Drying Time Input for Solid and Liquid Cycle

Intuitive Interface: Navigate with ease using the PLC (Programmable Logic Controller) paired with an HMI (Human Machine Interface) touch screen, providing user-friendly control at your fingertips.



Sterilization Progress Counter

Precision Temperature Management: Temperature-dependent control ensures accuracy and precision throughout the sterilization process, guaranteeing optimal results.




Sterilization Pause and Stop Function

Advanced Sterilization Technology: Benefit from the innovative S.F.P.P. (Steam Flush Pressure Pulse) control system, enhancing sterilization efficacy and performance.




Numerical Keypad for Input Value

Seamless Data Documentation: Keep track of sterilization cycles effortlessly with the built-in printer temperature data logger, providing comprehensive documentation for compliance and quality assurance.



Instructional Window and Warning Display

Integrated Steam Generation: Enjoy the convenience of an integrated steam generator, ensuring a consistent and reliable steam supply for uninterrupted operation.



Industrial Grade

These autoclaves are meticulously crafted for rapid and thorough steam sterilization of heat and moisture-resistant materials. Boasting advanced PLC and HMI touch screen control, the EMC4 series offers customizable pre-vacuum, gravity, and liquid test cycles.

TALK TO US:

8749 5915 | 8711 7292 | 8742 5090

esphar_manila@yahoo.com

1425 Rizal Avenue, Sta. Cruz, Manila

Customizable Volume Capacity:
Tailored to Customer Preferences within Manufacturing Guidelines

