

17*17 / 14*17 inch Csl+a-Si+TFT




MX-FPD3543H
Wired FPD (Front)



MX-FPD3543H
Wired FPD (Back)



MX-FPD4343 SIZE



MX-FPD4343WL
Wireless FPD (Front)



MX-FPD4343WL
Wireless FPD (Back)



Computer and DR Software
(Optional)




Laptop and DR Software
(Optional)

Parameters

Parameter	MX-FPD3543H	MX-FPD3543WL	MX-FPD4343H	MX-FPD4343WL
Type	α-Si + Csl			
Image size(cm)	35*43cm 14*17 Inch	35*43cm 14*17 Inch	43*43cm 17*17 Inch	43*43cm 17*17 Inch
Pixel pitch(μm)	140μm			
Pixels Matrix	2560*3072	2560*3072	3072*3072	3072*3072
A/D(bit)	16bit			
spatial resolution	3.6 LP/mm	3.6 LP/mm	3.6 LP/mm	3.6 LP/mm
Water tightness	IP54	IP54	IP54	IP54
Battery stand by	No	7H	No	7H
Weight	3.0kg	3.0kg	3.7kg	3.7kg
Dimensions (cm)	38.3*46*1.5	38.3*46*1.5	46*46*1.5	46*46*1.5
Wireless	No	Yes	No	Yes

Excellent Images



【pelvis for adult,85kg】
【80KV 32mAs】



【Foot for adult,70kg】
【60KV 10mAs】



【hand for adult,80kg】
【50KV 6.3mAs】