

Infant Incubator MCG2005



Free-step mattress tilting adjust



Phototherapy unit(option)



Removable for cleaning



Phototherapy unit(option)

Product feature

1. Air temp. and skin temp. servo-controlled by computer
2. $>37^{\circ}\text{C}$ temperature setting
3. Various self-check alarms
4. Free-step mattress tilting adjustment
5. Removable humidity reservoir, easy to clean
6. With the transfusion shelf and tray
7. Added with guardrail
8. RS-232 connector
9. Phototherapy(option)

Specification

Electrical requirement	$\sim 220\text{V } 50\text{Hz}$ or $\sim 230\text{V } 60\text{Hz}$
Power consumption	$\leq 450\text{VA}$
Air mode temp. range	$25.0^{\circ}\text{C} \sim 37.0^{\circ}\text{C}$
Skin mode temp. range	$34.0^{\circ}\text{C} \sim 37.0^{\circ}\text{C}$
Air over temp. range	$< 38.0^{\circ}\text{C}$
Air mode high temp. range	$37.1^{\circ}\text{C} \sim 39.0^{\circ}\text{C}$
Skin mode high temp. range	$37.1^{\circ}\text{C} \sim 38.0^{\circ}\text{C}$
High temp. mode alarm	$< 40.0^{\circ}\text{C}$
Air mode alarm	$\pm 3.0^{\circ}\text{C}$
Skin mode alarm	$\pm 1.0^{\circ}\text{C}$
Skin sensor accuracy	$\leq \pm 0.3^{\circ}\text{C}$
Air control accuracy	$\leq 0.5^{\circ}\text{C}$
Mattress temp. uniformity	$\leq 0.8^{\circ}\text{C}$