




MCG2013 Infant Radiant Warmer

- Skin temp. and manual control by micro-computer
- Over temperature alarm
- Available adjustment of radiator and mattress tilting
- With one illumination light
- APGAR timer
- RS-232 connector



Phototherapy unit(option)



Low-pressure suction



Radiation head that can be adjusted in horizontal angulation



Free-step mattress tilting adjust



X-ray tray(option)

Specifications Infant Incubator

Model	MCG2013
Electrical requirement	~ 220V 50Hz or ~ 230V 60Hz
Power consumption	≤ 700VA
Air mode temp. range	25.0℃ ~ 37.0℃
Over temp. alarm	< 38.0℃
High mode temp. range	37.1℃ ~ 38.0℃
High mode temp. alarm	< 39.0℃
Temp. deviation alarm	± 1.0℃
Skin sensor accuracy	≤ ± 0.3℃
Mattress temp. uniformity	≤ 2.0℃
Heating rate	≤ 30min
APGAR timing	1min、5min、10min

Model	MCG2013
screen	LED
skin temperature sensor	●
suction	option
O ₂ bottle	option
phototherapy	-
X-ray tray	●