

MCS1398 Ceiling Double Arm LED OR Light



Detail Parameter:

Illuminance (Lux)	40,000-160,000	40,000-160,000
Lamp bulb quantity	61 pcs	61 pcs
Brand of bulb	OSRAM(German)	
Bulb life	>50,000 hours	
Color Temp(K)	3700-5000	
Color rendering index(Ra)	85-98	
Light beam depth	120cm//47.2inch	
Diameter of spot	16-28cm//6.3-11inch	
Light adjust range	1%-100%	
Temperature (Operator head)	<15°C	
Input power	AC100-240V,50/60HZ 2.1A	
Output	24V 4.5A	
Best install height	2.7-3.1m	

Packing List:

NO	Item	Quantity
1	LED 5 head	1 UNIT
2	LED 5 head	1 UNIT
3	Rotating arm + fixing base	1 SET
4	Balance Arm	2 SETS
5	Switching power supply	2 UNITS
6	Sterilizer handle	4 PIECES
7	Allen Wrench	1 SET
8	Mounting fixing bolt	1 SET
9	Shroud Base	1 SET
10	Large shield	1 SET
11	Manual book	1 UNIT

MCS0123 Ceilling Double Arm LED OR Light



Detail Parameter:

Illuminancy(adjustable)	40,000-180,000Lux	30,000-160,000Lux
Lamp Bulb Quantity	61	39
Bulb Type	German OSRAM LED bulb	German OSRAM LED bulb
Bulb Life (h)	≥80, 000h	≥80, 000h
Color temp (Kelvin)	3700K ~5000K	3700K ~5000K
Color rendering index(Ra)	85~98	85~98
Light Beam Depth	50-180cm	50-180cm
Diameter of Spot (mm)	160-280mm	160-280mm
Light Adjust Range	1%-100%	1%-100%
Consumed power	65	40
Temp Rise (doctor head)	<1°C	<1°C
Input voltage	AC100-240V 50/60HZ	AC100-240V 50/60HZ
Best Install Height(m)	2.7-3.1m	2.7-3.1m

Packing List:

NO	Item	Quantity
1	LED 5 head	1 UNIT
2	LED 3 head	1 UNIT
3	Rotating arm + fixing base	1 SET
4	Balance Arm	2 SETS
5	Switching power supply	2 UNITS
6	Sterilizer handle	4 PIECES
7	Allen Wrench	1 SET
8	Mounting fixing bolt	1 SET
9	Shroud Base	1 SET
10	Large shield	1 SET
11	Manual book	1 UNIT

MCS0500 Ceilling Single Arm LED OR Light



Illuminancy(adjustable)	40,000-180,000Lux
Lamp Bulb Quantity	61
Bulb Type	German OSRAM LED bulb
Bulb Life (h)	≥80, 000h
Color temp (Kelvin)	3700K ~5000K
Color rendering index(Ra)	85~98
Light Beam Depth	50-180cm
Diameter of Spot (mm)	160-280mm
Light Adjust Range	1%-100%
Consumed power	65
Temp Rise (doctor head)	<1°C
Input voltage	AC100-240V 50/60HZ
Best Install Height(m)	2.7-3.1m

Packing List:

NO	Item	Quantity
1	LED 5 head	1 UNIT
3	Rotating arm + fixing base	1 SET
4	Balance Arm	1 SETS
5	Switching power supply	1 UNITS
6	Sterilizer handle	2 PIECES
7	Allen Wrench	1 SET
8	Mounting fixing bolt	1 SET
9	Shroud Base	1 SET
10	Large shield	1 SET
11	Manual book	1 UNIT

MCS0500 Mobile LED Operation Lamp



Detail Parameter:

Illuminance (Lux)	40,000-160,000
Lamp bulb quantity	61 pcs
Brand of bulb	OSRAM(German)
Bulb life	>50,000 hours
Color Temp(K)	3700-5000
Color rendering index(Ra)	85-98
Lightbeamdepeth	120cm//47.2inch
Diameter of spot	16-28cm//6.3-11inch
Light adjust range	1%-100%
Temp rise(Operator head)	<1.5°C
Input power	AC100-240V,50/60HZ 2.1A
Output	24V 4.5A
Best install height	2.7-3.1m

Packing List:

NO	Item	Quantity
1	LED 5 head	1 UNIT
2	Rotating arm + fixing base	1 SET
3	Balance Arm	1 SETS
4	Switching power supply	1 UNITS
5	Sterilizer handle	2 PIECES
6	Allen Wrench	1 SET
7	Mounting fixing bolt	1 SET
8	Shroud Base	1 SET
9	Large shield	1 SET
10	Manual book	1 UNIT

