

Operating Table

Model: MCS0646



Instructions:

MCS0646 electric operating table is used for general surgery, heart and kidney, orthopedics, neurosurgery, gynecology, urology and other surgical. It's driven by low-voltage DC motor with advanced backup power, safe and reliable. The leg board can open and is removable. The top board can move back and forth by 300mm, using with the C-arm.

The base cover is made of stainless steel, slushing and easy to clean.

Main technical parameters:

Length	2020mm
Width	500mm
Min height	660mm
Max height	960mm
Table translation	300mm